

GSK3B / GSK3 Beta Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C219605 - LSBio	
CatalogID:	LS-C219605
Target:	glycogen synthase kinase 3 beta (GSK3B)
Synonyms:	GSK3B Antibody, GSK-3 beta Antibody, GSK3 beta Antibody, GSK3beta Antibody, GSK3beta isoform Antibody, TPKI Antibody, GSK 3beta Antibody, Tau protein kinase I Antibody, ZGSK-3beta Antibody
Family / Subfamily:	Protein Kinase / GSK-3
Host	GSK3B antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	GSK3B / GSK3 Beta antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GSK3B / GSK3 Beta antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human GSK3B.
Specificity:	Human GSK3B / GSK3 Beta
Epitope:	N-Terminus
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	