

HSPB1 / HSP27 Rabbit anti-Human Polyclonal (Ser78) (FITC) Antibody - LS-C219679 - LSBio

<b>CatalogID:</b>	LS-C219679
<b>Target:</b>	heat shock 27kDa protein 1 (HSPB1)
<b>Synonyms:</b>	HSPB1 Antibody, 28 kDa heat shock protein Antibody, Heat shock 27 kDa protein Antibody, Heat shock 27kD protein 1 Antibody, Heat shock 27kDa protein 1 Antibody, HSP27 Antibody, HSP28 Antibody, HSP 27 Antibody, Hsp25 Antibody, SRP27 Antibody, Stress-responsive protein 27 Antibody, CMT2F Antibody, Heat shock 27 kd protein Antibody, Heat shock protein beta-1 Antibody, HMN2B Antibody, HS.76067 Antibody
<b>Host</b>	HSPB1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	HSPB1 / HSP27 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	HSPB1 / HSP27 antibody was raised against kLH-conjugated synthetic non-phosphopeptide around Ser78 of human HSPB1.
<b>Specificity:</b>	Human HSPB1 / HSP27
<b>Epitope:</b>	Ser78
<b>Reactivity:</b>	Human
<b>Purification:</b>	Immunoaffinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences