

IL-1Rrp2 / IL1RL2 Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C219776 - LSBio

CatalogID:	LS-C219776
Target:	interleukin 1 receptor-like 2 (IL1RL2)
Synonyms:	IL1RL2 Antibody, IL-1Rrp2 Antibody, IL1R-rp2 Antibody, IL1RRP2 Antibody, Interleukin 1 receptor-like 2 Antibody, Interleukin-1 receptor-like 2 Antibody, IL-36R Antibody
Family / Subfamily:	Interleukin
Host	IL1RL2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IL-1Rrp2 / IL1RL2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	IL-1Rrp2 / IL1RL2 antibody was raised against kLH-conjugated synthetic peptide from the C-terminal(E487) region of human IL1R.
Specificity:	Human IL-1Rrp2 / IL1RL2
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences