

IL2 Rabbit anti-Human Polyclonal (Internal) (FITC) Antibody - LS-C219803 - LSBio

CatalogID:	LS-C219803
Target:	interleukin 2 (IL2)
Synonyms:	IL2 Antibody, Aldesleukin Antibody, Interleukin 2 Antibody, Lymphokine Antibody, IL-2 Antibody, Interleukin-2 Antibody, TCGF Antibody, T cell growth factor Antibody, T-cell growth factor Antibody
Family / Subfamily:	Interleukin
Host	IL2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IL2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	IL2 antibody was raised against kLH-conjugated synthetic peptide from internal region of human IL2.
Specificity:	Human IL2 / IL-2
Epitope:	Internal
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences