

NME2 / NM23 Rabbit anti-Human Polyclonal (N-Terminus) (AP) Antibody - LS-C220189 - LSBio	
CatalogID:	LS-C220189
Target:	NME/NM23 nucleoside diphosphate kinase 2 (NME2)
Synonyms:	NME2 Antibody, Histidine protein kinase NDKB Antibody, NDKB Antibody, NM23-H2 Antibody, PUF Antibody, NDP kinase B Antibody, NDPK-B Antibody, C-myc transcription factor Antibody, NDK B Antibody, NDPKB Antibody, NM23 Antibody, NM23B Antibody
Host	NME2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Alkaline Phosphatase (AP)
Immunogen Species:	NME2 / NM23 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	NME2 / NM23 antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human NME2.
Specificity:	Human NME2 / NM23.
Epitope:	N-Terminus
Reactivity:	Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	