

Osteocalcin Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C220217 - LSBio	
CatalogID:	LS-C220217
Target:	bone gamma-carboxyglutamate (gla) protein (BGLAP)
Synonyms:	BGLAP Antibody, Bone Gla protein Antibody, BGP Antibody, OC Antibody, Osteocalcin Antibody
Host	BGLAP antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Osteocalcin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	Osteocalcin antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human osteocalcin.
Specificity:	Human Osteocalcin.
Epitope:	N-Terminus
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	