

AGER / RAGE Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C220517 - LSBio	
CatalogID:	LS-C220517
Target:	advanced glycosylation end product-specific receptor (AGER)
Synonyms:	AGER Antibody, RAGE Antibody, RAGE isoform NtRAGE-delta Antibody, RAGE isoform sRAGE-delta Antibody
Host	AGER antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	AGER / RAGE antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AGER / RAGE antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human AGER.
Specificity:	Human AGER / RAGE.
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	