

RGS19 Rabbit anti-Human Polyclonal (pSer24) (PE) Antibody - LS-C220555 - LSBio

CatalogID:	LS-C220555
Target:	regulator of G-protein signaling 19 (RGS19)
Synonyms:	RGS19 Antibody, GNAI3IP Antibody, G alpha interacting protein Antibody, G-alpha-interacting protein Antibody, RGSGAIP Antibody, GAIP Antibody
Family / Subfamily:	Regulator of G-protein signalling / not assigned-Regulator of G-protein signalling
Host	RGS19 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	RGS19 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	RGS19 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser24 of human RGS19.
Specificity:	Human RGS19
Epitope:	pSer24
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences