

SMAD3 Rabbit anti-Human Polyclonal (pSer208) (PE) Antibody - LS-C220703 - LSBio

CatalogID:	LS-C220703
Target:	SMAD family member 3 (SMAD3)
Synonyms:	SMAD3 Antibody, HSPC193 Antibody, HMAD-3 Antibody, HsT17436 Antibody, LDS1C Antibody, LDS3 Antibody, Mad homolog JV15-2 Antibody, MAD3 Antibody, MADH3 Antibody, Mothers against DPP homolog 3 Antibody, JV15-2 Antibody, SMA- and MAD-related protein 3 Antibody, SMAD family member 3 Antibody, MAD homolog 3 Antibody, HSMAD3 Antibody, Mad protein homolog Antibody, SMAD 3 Antibody
Host	SMAD3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	SMAD3 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	SMAD3 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser208 of human SMAD3.
Specificity:	Human SMAD3.
Epitope:	pSer208
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ICC, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences