

Vimentin Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C221237 - LSBio

CatalogID:	LS-C221237
Target:	vimentin
Synonyms:	VIM Antibody, Vimentin Antibody
Host	VIM antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Vimentin antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	Vimentin antibody was raised against recombinant human vimentin.
Specificity:	Human Vimentin.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences