

ABI2B / ABI2 Rabbit anti-Human Polyclonal (aa1-30) (FITC) Antibody - LS-C221320 - LSBio

CatalogID:	LS-C221320
Target:	abl-interactor 2 (ABI2)
Synonyms:	ABI2 Antibody, Abelson interactor 2 Antibody, Abl interactor 2 Antibody, Abl-binding protein 3 Antibody, ArgBP1 Antibody, ArgBPIB Antibody, Arg-binding protein 1 Antibody, Argbpib protein Antibody, ABI-2 Antibody, ABI2B Antibody, Abl binding protein 3 Antibody, Abl-interactor 2 Antibody, Abl-interactor protein 2b Antibody, AbIBP3 Antibody, AIP-1 Antibody, SSH3BP2 Antibody, ArgBPIA Antibody
Host	ABI2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ABI2B / ABI2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ABI2B / ABI2 antibody was raised against kLH-conjugated synthetic peptide (aa1-30) from the N-terminal region of human ABI2.
Specificity:	Human ABI2B / ABI2
Epitope:	aa1-30
Reactivity:	Human, Mouse
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences