

ARFGAP1 Rabbit anti-Human Polyclonal (aa377-406) (FITC) Antibody - LS-C221461 - LSBio	
CatalogID:	LS-C221461
Target:	ADP-ribosylation factor GTPase activating protein 1 (ARFGAP1)
Synonyms:	ARFGAP1 Antibody, ARF GAP 1 Antibody, ARF1GAP Antibody, ARF1 GAP Antibody, BA261N11.3 Antibody, GAP protein Antibody, HRIHFB2281 Antibody
Host	ARFGAP1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ARFGAP1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ARFGAP1 antibody was raised against kLH-conjugated synthetic peptide (aa377-406) from the C-terminus of human ARFGAP1.
Specificity:	Human ARFGAP1.
Epitope:	aa377-406
Reactivity:	Human, Rabbit
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	