

ATL3 Rabbit anti-Human Polyclonal (aa230-260) (FITC) Antibody - LS-C221503 - LSBio

CatalogID:	LS-C221503
Target:	atlastin GTPase 3 (ATL3)
Synonyms:	ATL3 Antibody, Atlastin GTPase 3 Antibody, Atlastin-3 Antibody, DKFZP564J0863 Antibody, Atlastin 3 Antibody
Host	ATL3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ATL3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ATL3 antibody was raised against kLH-conjugated synthetic peptide between amino acids 230-260 from human ATL3.
Specificity:	Human ATL3
Epitope:	aa230-260
Reactivity:	Human, Mouse, Hamster
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences