

BID Rabbit anti-Human Polyclonal (aa78-108) (PE) Antibody - LS-C221535 - LSBio

CatalogID:	LS-C221535
Target:	BH3 interacting domain death agonist (BID)
Synonyms:	BID Antibody, Apoptic death agonist Antibody, BID isoform L(2) Antibody, BID isoform Si6 Antibody, Bh3 Antibody, FP497 Antibody, p22 BID Antibody, BID isoform ES(1b) Antibody, Desmocollin type 4 Antibody, Human BID coding sequence Antibody
Family / Subfamily:	Apoptosis
Host	BID antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	BID antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	BID antibody was raised against KLH-conjugated synthetic peptide (aa78-108) from the BH3 domain of human Bid.
Specificity:	Human BID.
Epitope:	aa78-108
Reactivity:	Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences