

BMI1 / PCGF4 Rabbit anti-Human Polyclonal (aa128-158) (FITC) Antibody - LS-C221539 - LSBio	
CatalogID:	LS-C221539
Target:	BMI1 polycomb ring finger oncogene
Synonyms:	BMI1 Antibody, BMI-1 Antibody, Flvi-2/bmi-1 Antibody, FLVI2/BMI1 Antibody, PCGF4 Antibody, Polycomb group protein Bmi1 Antibody, Polycomb complex protein BMI-1 Antibody, Dna-binding protein bmi-1 Antibody, Polycomb group ring finger 4 Antibody, RING finger protein 51 Antibody, RNF51 Antibody
Host	BMI1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	BMI1 / PCGF4 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	BMI1 / PCGF4 antibody was raised against kLH-conjugated synthetic peptide (aa128-158) of human BMI.
Specificity:	Human BMI1 / BMI-1
Epitope:	aa128-158
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014 © 2014 LifeSpan BioSciences	