

CD105 Rabbit anti-Human Polyclonal (aa387-417) (FITC) Antibody - LS-C221690 - LSBio	
CatalogID:	LS-C221690
Target:	endoglin (ENG)
Synonyms:	ENG Antibody, CD105 Antibody, CD105 antigen Antibody, END Antibody, HHT1 Antibody, ORW Antibody, ORW1 Antibody, Osler-Rendu-Weber syndrome 1 Antibody, Endoglin Antibody
Host	ENG antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD105 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CD105 antibody was raised against kLH-conjugated synthetic peptide (aa387-417) from internal region of human CD105.
Specificity:	Human CD105.
Epitope:	aa387-417
Reactivity:	Human, Mouse
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	