

CD69 Rabbit anti-Human Polyclonal (aa126-155) (FITC) Antibody - LS-C221708 - LSBio

CatalogID:	LS-C221708
Target:	CD69 molecule
Synonyms:	CD69 Antibody, Activation inducer molecule Antibody, AIM Antibody, CD69 antigen Antibody, CLEC2C Antibody, EA1 Antibody, Early activation antigen CD69 Antibody, BL-AC/P26 Antibody, CD69 molecule Antibody, gp32/28 Antibody, MLR-3 Antibody
Host	CD69 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD69 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CD69 antibody was raised against kLH-conjugated synthetic peptide (aa126-155) of human CD69.
Specificity:	Human CD69
Epitope:	aa126-155
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences