

CLDN7 / Claudin 7 Rabbit anti-Human Polyclonal (aa179-209) (Biotin) Antibody - LS-C221806 - LSBio	
<b>CatalogID:</b>	LS-C221806
<b>Target:</b>	claudin 7 (CLDN7)
<b>Synonyms:</b>	CLDN7 Antibody, Claudin-7 Antibody, Claudin-1 Antibody, CLDN-7 Antibody, Hs.84359 Antibody, CEPTRL2 Antibody, Claudin 7 Antibody, CPETRL2 Antibody
<b>Family / Subfamily:</b>	Claudin / not assigned-Claudin
<b>Host</b>	CLDN7 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Biotin
<b>Immunogen Species:</b>	CLDN7 / Claudin 7 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	CLDN7 / Claudin 7 antibody was raised against kLH-conjugated synthetic peptide (aa179-209) from the C-terminus of human CLDN7.
<b>Specificity:</b>	Human CLDN7 / Claudin 7
<b>Epitope:</b>	aa179-209
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Uses:</b>	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b> <b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b> Created on 9/28/2014 © 2014 LifeSpan BioSciences	