

CLDN7 / Claudin 7 Rabbit anti-Human Polyclonal (aa179-209) (PE) Antibody - LS-C221813 - LSBio	
CatalogID:	LS-C221813
Target:	claudin 7 (CLDN7)
Synonyms:	CLDN7 Antibody, Claudin-7 Antibody, Claudin-1 Antibody, CLDN-7 Antibody, Hs.84359 Antibody, CEPTRL2 Antibody, Claudin 7 Antibody, CPETRL2 Antibody
Family / Subfamily:	Claudin / not assigned-Claudin
Host	CLDN7 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CLDN7 / Claudin 7 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CLDN7 / Claudin 7 antibody was raised against kLH-conjugated synthetic peptide (aa179-209) from the C-terminus of human CLDN7.
Specificity:	Human CLDN7 / Claudin 7
Epitope:	aa179-209
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	