

C1orf107 / DEF Rabbit anti-Human Polyclonal (aa9-38) (PE) Antibody - LS-C222060 - LSBio

<b>CatalogID:</b>	LS-C222060
<b>Target:</b>	digestive organ expansion factor homolog (zebrafish) (DIEXF)
<b>Synonyms:</b>	DIEXF Antibody, C1orf107 Antibody, DJ434O14.5 Antibody, UTP25 Antibody, DEF Antibody
<b>Host</b>	DIEXF antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	C1orf107 / DEF antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	C1orf107 / DEF antibody was raised against kLH-conjugated synthetic peptide (aa9-38) from the N-terminal region of human DEF.
<b>Specificity:</b>	Human C1orf107 / DEF.
<b>Epitope:</b>	aa9-38
<b>Reactivity:</b>	Human, Mouse, Hamster
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Uses:</b>	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences