

DEK Rabbit anti-Human Polyclonal (aa343-372) (FITC) Antibody - LS-C222064 - LSBio

CatalogID:	LS-C222064
Target:	DEK oncogene
Synonyms:	DEK Antibody, D6S231E Antibody, DEK oncogene (DNA binding) Antibody, Protein DEK Antibody, DEK oncogene Antibody
Host	DEK antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	DEK antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	DEK antibody was raised against kLH-conjugated synthetic peptide (aa343-372) from the C-terminus of human DEK.
Specificity:	Human DEK
Epitope:	aa343-372
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences