

DGKZ Rabbit anti-Human Polyclonal (aa997-1026) (PE) Antibody - LS-C222084 - LSBio

CatalogID:	LS-C222084
Target:	diacylglycerol kinase, zeta (DGKZ)
Synonyms:	DGKZ Antibody, DAG kinase zeta Antibody, DAGK6 Antibody, Diglyceride kinase zeta Antibody, Diacylglycerol kinase zeta Antibody, DGK-ZETA Antibody, Diacylglycerol kinase, zeta Antibody, HDGKzeta Antibody, DAGK5 Antibody
Family / Subfamily:	Non-protein kinase / Diacylglycerol kinase
Host	DGKZ antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	DGKZ antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	DGKZ antibody was raised against kLH-conjugated synthetic peptide (aa997-1026) from the C-terminus of human DGKZ.
Specificity:	Human DGKZ.
Epitope:	aa997-1026
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences