

DNMT / DNMT1 Rabbit anti-Human Polyclonal (pSer714) (FITC) Antibody - LS-C222180 - LSBio	
<b>CatalogID:</b>	LS-C222180
<b>Target:</b>	DNA (cytosine-5-)-methyltransferase 1 (DNMT1)
<b>Synonyms:</b>	DNMT1 Antibody, AIM Antibody, CXXC finger protein 9 Antibody, CXXC9 Antibody, DNA methyltransferase 1 Antibody, DNMT Antibody, DNA methyltransferase Hsa1 Antibody, DNA MTase Hsa1 Antibody, HSN1E Antibody, M.Hsa1 Antibody, MCMT Antibody
<b>Family / Subfamily:</b>	C5-methyltransferase
<b>Host</b>	DNMT1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	DNMT / DNMT1 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	DNMT / DNMT1 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser714 of human Dnmt1.
<b>Specificity:</b>	Human DNMT / DNMT1
<b>Epitope:</b>	pSer714
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Uses:</b>	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
<b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b>	
Created on 9/28/2014 © 2014 LifeSpan BioSciences	