

DNMT / DNMT1 Rabbit anti-Human Polyclonal (pSer714) (PE) Antibody - LS-C222182 - LSBio

CatalogID:	LS-C222182
Target:	DNA (cytosine-5-)-methyltransferase 1 (DNMT1)
Synonyms:	DNMT1 Antibody, AIM Antibody, CXXC finger protein 9 Antibody, CXXC9 Antibody, DNA methyltransferase 1 Antibody, DNMT Antibody, DNA methyltransferase Hsa1 Antibody, DNA MTase Hsa1 Antibody, HSN1E Antibody, M.Hsa1 Antibody, MCMT Antibody
Family / Subfamily:	C5-methyltransferase
Host	DNMT1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	DNMT / DNMT1 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	DNMT / DNMT1 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser714 of human Dnmt1.
Specificity:	Human DNMT / DNMT1
Epitope:	pSer714
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences