

DNMT / DNMT1 Rabbit anti-Human Polyclonal (pSer1105) (FITC) Antibody - LS-C222185 - LSBio	
CatalogID:	LS-C222185
Target:	DNA (cytosine-5-)-methyltransferase 1 (DNMT1)
Synonyms:	DNMT1 Antibody, AIM Antibody, CXXC finger protein 9 Antibody, CXXC9 Antibody, DNA methyltransferase 1 Antibody, DNMT Antibody, DNA methyltransferase Hsa1 Antibody, DNA MTase Hsa1 Antibody, HSN1E Antibody, M.Hsa1 Antibody, MCMT Antibody
Family / Subfamily:	C5-methyltransferase
Host	DNMT1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	DNMT / DNMT1 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	DNMT / DNMT1 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser1105 of human Dnmt1.
Specificity:	Human DNMT / DNMT1
Epitope:	pSer1105
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	