

DPP4 / CD26 Rabbit anti-Human Polyclonal (aa46-74) (FITC) Antibody - LS-C222229 - LSBio	
CatalogID:	LS-C222229
Target:	dipeptidyl-peptidase 4 (DPP4)
Synonyms:	DPP4 Antibody, ADCP-2 Antibody, ADABP Antibody, Dipeptidyl peptidase IV Antibody, Dipeptidyl-peptidase 4 Antibody, DPPIV Antibody, ADCP2 Antibody, T-cell activation antigen CD26 Antibody, CD26 Antibody, CD26 antigen Antibody, Dipeptidyl peptidase 4 Antibody, Dipeptidylpeptidase 4 Antibody, DPP IV Antibody, TP103 Antibody
Family / Subfamily:	Protease / Serine S9B
Host	DPP4 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	DPP4 / CD26 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	DPP4 / CD26 antibody was raised against kLH-conjugated synthetic peptide (aa46-74) from the N-terminal region of human DPP4.
Specificity:	Human DPP4 / CD26.
Epitope:	aa46-74
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014 © 2014 LifeSpan BioSciences	