

ETFA Rabbit anti-Human Polyclonal (aa283-311) (PE) Antibody - LS-C222477 - LSBio

CatalogID:	LS-C222477
Target:	electron-transfer-flavoprotein, alpha polypeptide (ETFA)
Synonyms:	ETFA Antibody, Alpha-ETF Antibody, GA2 Antibody, EMA Antibody, Glutaric aciduria II Antibody
Host	ETFA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ETFA antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ETFA antibody was raised against kLH-conjugated synthetic peptide (aa283-311) from the C-terminus of human ETFA.
Specificity:	Human ETFA
Epitope:	aa283-311
Reactivity:	Human, Mouse
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences