

CLTA / LCA Chicken anti-Human Polyclonal (aa1-211) Antibody - LS-C22250 - LSBio

CatalogID:	LS-C22250
Target:	clathrin, light chain A (CLTA)
Synonyms:	CLTA Antibody, Clathrin light chain A Antibody, Clathrin, light chain A Antibody, LCA Antibody
Host	CLTA antibody was produced in Chicken
Clonality:	Polyclonal
Isotype:	IgY
Immunogen Species:	CLTA / LCA antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CLTA / LCA antibody was raised against synthetic peptide corresponding to aa1-211.
Specificity:	Recognizes human Clathrin, Light Chain, isoform a.
Epitope:	aa1-211
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.02% sodium azide.
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Usage Summary:	Suitable for use in ELISA and Western Blot. Western Blot: 1:1000-1:2000.
Uses:	Western blot (1:1000 - 1:2000), ELISA (Optimal dilution to be determined by the researcher)
Size:	50 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences