

FCGR3A / CD16A Rabbit anti-Human Polyclonal (aa225-254) (FITC) Antibody - LS-C222543 - LSBio

CatalogID:	LS-C222543
Target:	Fc fragment of IgG, low affinity IIIa, receptor (CD16a) (FCGR3A)
Synonyms:	FCGR3A Antibody, CD16A Antibody, CD16a antigen Antibody, CD16 Antibody, Fc gammaRIIIA Antibody, FCR-10 Antibody, FCRIII Antibody, Fc gamma receptor III-A Antibody, Fc-gamma RIII Antibody, Fc-gamma RIII-alpha Antibody, FCG3 Antibody, IGFR3 Antibody, Fc- gamma riii Antibody, Fc-gamma receptor IIIb (CD16) Antibody, Fc-gamma RIIIa Antibody, FCGR3 Antibody, FCRIII Antibody, Neutrophil-specific antigen NA Antibody, Igg fc receptor iii-1 Antibody, IgG Fc receptor III-2 Antibody, Igg fc receptors iii-1 Antibody, FCGRIII Antibody
Host	FCGR3A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	FCGR3A / CD16A antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	FCGR3A / CD16A antibody was raised against kLH-conjugated synthetic peptide (aa225-254) from the C-terminus of human FCGR3A.
Specificity:	Human FCGR3A / CD16A
Epitope:	aa225-254
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences