

| FKRP Rabbit anti-Human Polyclonal (aa236-265) (FITC) Antibody - LS-C222561 - LSBio  |   |
|---|---|
| <b>CatalogID:</b>   | LS-C222561  |
| <b>Target:</b>  | fukutin related protein (FKRP)  |
| <b>Synonyms:</b>  | FKRP Antibody, Fukutin related protein Antibody, Fukutin-related protein Antibody, MDC1C Antibody, MDDGB5 Antibody, MDDGC5 Antibody, LGMD2I Antibody, MDDGA5 Antibody |
| <b>Host</b>   | FKRP antibody was produced in Rabbit  |
| <b>Clonality:</b>   | Polyclonal  |
| <b>Isotype:</b>   | IgG   |
| <b>Conjugations:</b>  | Fluorescein (FITC)  |
| <b>Immunogen Species:</b>   | FKRP antibody was raised against Human  |
| <b>Antigen Type:</b>  | Synthetic peptide   |
| <b>Immunogen:</b>   | FKRP antibody was raised against kLH-conjugated synthetic peptide (aa236-265) of human FKRP.  |
| <b>Specificity:</b>   | Human FKRP  |
| <b>Epitope:</b>   | aa236-265   |
| <b>Reactivity:</b>  | Human   |
| <b>Purification:</b>  | Protein A purified  |
| <b>Presentation:</b>  | PBS, pH 7.2, 0.09% sodium azide.  |
| <b>Recommended Storage:</b>   | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.                           |
| <b>Uses:</b>  | Western blot, ELISA (Optimal dilution to be determined by the researcher)   |
| <b>Size:</b>  | 200 µl  |
| <b>Requested From:</b>  | Japan   |
| <b>Laboratory Reagent For In Vitro Research Use Only</b>                            |   |
| <b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b> |   |
| Created on 9/28/2014<br>© 2014 LifeSpan BioSciences                                 |   |