

AAGAB Rabbit anti-Human Polyclonal (aa293-322) (PE) Antibody - LS-C222569 - LSBio

CatalogID:	LS-C222569
Target:	alpha- and gamma-adaptin binding protein (AAGAB)
Synonyms:	AAGAB Antibody, p34 Antibody, PPKP1A Antibody, PPKP1 Antibody, FLJ11506 Antibody, KPPP1 Antibody
Host	AAGAB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	AAGAB antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AAGAB antibody was raised against kLH-conjugated synthetic peptide (aa293-322) from the C-terminus of human FLJ11506.
Specificity:	Human AAGAB.
Epitope:	aa293-322
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences