

Gephyrin Rabbit anti-Human Polyclonal (aa444-474) (FITC) Antibody - LS-C222648 - LSBio

CatalogID:	LS-C222648
Target:	gephyrin (GPHN)
Synonyms:	GPHN Antibody, Gephyrin Antibody, KIAA1385 Antibody, GEPH Antibody, GPH Antibody, GPHRYN Antibody, HKPX1 Antibody
Host	GPHN antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Gephyrin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	Gephyrin antibody was raised against kLH-conjugated synthetic peptide (aa444-474) of human GPHN.
Specificity:	Human Gephyrin
Epitope:	aa444-474
Reactivity:	Human, Mouse
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences