

HIAT1 Rabbit anti-Human Polyclonal (aa460-490) (PE) Antibody - LS-C222734 - LSBio

CatalogID:	LS-C222734
Target:	hippocampus abundant transcript 1 (HIAT1)
Synonyms:	HIAT1 Antibody, DKFZP564L0864 Antibody
Family / Subfamily:	Transporter / not assigned-Transporter
Host	HIAT1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	HIAT1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	HIAT1 antibody was raised against kLH-conjugated synthetic peptide between amino acids 460-490 from the C-terminus of human HIAT1.
Specificity:	Human HIAT1
Epitope:	aa460-490
Reactivity:	Human, Mouse, Hamster
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences