

Complement C1QA Sheep anti-Human Polyclonal (FITC) Antibody - LS-C22274 - LSBio

CatalogID:	LS-C22274
Target:	complement component 1, q subcomponent, A chain (C1QA)
Synonyms:	C1QA Antibody, Complement C1QA Antibody
Host	C1QA antibody was produced in Sheep
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Complement C1QA antibody was raised against Human
Antigen Type:	Purified protein
Immunogen:	Complement C1QA antibody was raised against human C1q, purified from plasma.
Specificity:	Recognizes human C1q, a sub-component of the C1 complement component.
Reactivity:	Human
Purification:	Ion exchange chromatography
Presentation:	PBS, pH 7.2, 0.1% sodium azide.
Usage Summary:	Immunohistochemistry: (frozen)- (direct) 1:25-1:50.
Uses:	IHC - Frozen (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences