

HSP90B1 / GP96 / GRP94 Rabbit anti-Human Polyclonal (aa20-50) (PE) Antibody - LS-C222768 - LSBio	
CatalogID:	LS-C222768
Target:	heat shock protein 90kDa beta (Grp94), member 1 (HSP90B1)
Synonyms:	HSP90B1 Antibody, ECGP Antibody, Endoplasmic Antibody, gp96 homolog Antibody, gp96 Antibody, GRP-94 Antibody, GRP94 Antibody, Tumor rejection antigen 1 Antibody, TRA1 Antibody
Host	HSP90B1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	HSP90B1 / GP96 / GRP94 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	HSP90B1 / GP96 / GRP94 antibody was raised against kLH-conjugated synthetic peptide between amino acids 20-50 from the N-terminal region of human HSP90B1.
Specificity:	Human HSP90B1 / GP96 / GRP94
Epitope:	aa20-50
Reactivity:	Human, Mouse, Hamster
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	