

MAGOH Rabbit anti-Human Polyclonal (aa45-74) (PE) Antibody - LS-C223133 - LSBio

CatalogID:	LS-C223133
Target:	mago-nashi homolog, proliferation-associated (Drosophila) (MAGOH)
Synonyms:	MAGOH Antibody, MAGOHA Antibody, MAGOH1 Antibody, Mago nashi protein homolog Antibody, Protein mago nashi homolog Antibody
Host	MAGOH antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	MAGOH antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MAGOH antibody was raised against KLH-conjugated synthetic peptide (aa45-74) of human MAGOH.
Specificity:	Human MAGOH.
Epitope:	aa45-74
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences