

MAP1LC3A / LC3A Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C223140 - LSBio

CatalogID:	LS-C223140
Target:	microtubule-associated protein 1 light chain 3 alpha (MAP1LC3A)
Synonyms:	MAP1LC3A Antibody, ATG8E Antibody, LC3 Antibody, LC3A Antibody, MAP1BLC3 Antibody, MAP1A/1B light chain 3 A Antibody, MAP1A/MAP1B LC3 A Antibody, MAP1A/MAP1B light chain 3 A Antibody, MAP1ALC3 Antibody
Host	MAP1LC3A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	MAP1LC3A / LC3A antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	MAP1LC3A / LC3A antibody was raised against recombinant human LC3.
Specificity:	Human MAP1LC3A / LC3A.
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences