

NDRG2 Rabbit anti-Human Polyclonal (aa41-70) (AP) Antibody - LS-C223421 - LSBio

CatalogID:	LS-C223421
Target:	NDRG family member 2 (NDRG2)
Synonyms:	NDRG2 Antibody, KIAA1248 Antibody, N-myc downstream regulator 2 Antibody, NDRG family member 2 Antibody, SYLD Antibody, NDR1-related protein NDR2 Antibody, Protein Syld709613 Antibody, Syld709613 protein Antibody, Cytoplasmic protein Ndr1 Antibody, Protein NDRG2 Antibody
Host	NDRG2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Alkaline Phosphatase (AP)
Immunogen Species:	NDRG2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	NDRG2 antibody was raised against kLH-conjugated synthetic peptide (aa41-70) from the N-terminal region of human NDRG2.
Specificity:	Human NDRG2
Epitope:	aa41-70
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences