

NFKBIL1 Rabbit anti-Human Polyclonal (aa256-289) (FITC) Antibody - LS-C223458 - LSBio	
<b>CatalogID:</b>	LS-C223458
<b>Target:</b>	nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor-like 1 (NFKBIL1)
<b>Synonyms:</b>	NFKBIL1 Antibody, I-kappa-B-like protein Antibody, NFKBIL Antibody, IkappaBL Antibody, IKBL Antibody
<b>Host</b>	NFKBIL1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	NFKBIL1 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	NFKBIL1 antibody was raised against kLH-conjugated synthetic peptide (aa256-289) of human NFKBIL1.
<b>Specificity:</b>	Human NFKBIL1
<b>Epitope:</b>	aa256-289
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Uses:</b>	ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	