

PFKP Rabbit anti-Human Polyclonal (aa755-784) (PE) Antibody - LS-C223633 - LSBio

CatalogID:	LS-C223633
Target:	phosphofructokinase, platelet (PFKP)
Synonyms:	PFKP Antibody, 6-phosphofructokinase type C Antibody, 6-phosphofructokinase, type c Antibody, PFK-C Antibody, PFKF Antibody, Phosphofructokinase 1 Antibody, Phosphofructokinase, platelet Antibody, Phosphohexokinase Antibody
Host	PFKP antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	PFKP antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PFKP antibody was raised against kLH-conjugated synthetic peptide (aa755-784) from the C-terminus of human PFKP.
Specificity:	Human PFKP
Epitope:	aa755-784
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences