

SORBS1 / Ponsin Rabbit anti-Human Polyclonal (aa869-896) (PE) Antibody - LS-C224118 - LSBio	
CatalogID:	LS-C224118
Target:	sorbin and SH3 domain containing 1 (SORBS1)
Synonyms:	SORBS1 Antibody, C-Cbl associated protein Antibody, C-Cbl-associated protein Antibody, CAP Antibody, FLAF2 Antibody, KIAA0894 Antibody, KIAA1296 Antibody, R85FL Antibody, SH3 domain protein 5 Antibody, SH3-domain protein 5 (ponsin) Antibody, SH3P12 Antibody, SORB1 Antibody, Fas-ligand associated factor 2 Antibody, Ponsin Antibody, SH3D5 Antibody
Host	SORBS1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	SORBS1 / Ponsin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	SORBS1 / Ponsin antibody was raised against kLH-conjugated synthetic peptide (aa869-896) of human SORBS1.
Specificity:	Human SORBS1 / Ponsin.
Epitope:	aa869-896
Reactivity:	Human, Hamster
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014 © 2014 LifeSpan BioSciences	