

VASP Rabbit anti-Human Polyclonal (aa266-296) (PE) Antibody - LS-C224627 - LSBio

CatalogID:	LS-C224627
Target:	vasodilator-stimulated phosphoprotein (VASP)
Synonyms:	VASP Antibody
Family / Subfamily:	Ena-VASP
Host	VASP antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	VASP antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	VASP antibody was raised against kLH-conjugated synthetic peptide (aa266-296) from the C-terminus of human VASP.
Specificity:	Human VASP.
Epitope:	aa266-296
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences