

**ATM Rabbit anti-Human Polyclonal (pSer1981) (FITC) Antibody - LS-C224863 - LSBio**

<b>CatalogID:</b>	LS-C224863
<b>Target:</b>	ataxia telangiectasia mutated (ATM)
<b>Synonyms:</b>	ATM Antibody, A-T mutated Antibody, ATA Antibody, ATC Antibody, ATD Antibody, ATE Antibody, AT mutated Antibody, Ataxia telangiectasia Antibody, ATDC Antibody, Serine-protein kinase ATM Antibody, AT1 Antibody, Ataxia telangiectasia mutated Antibody, TEL1 Antibody, TELO1 Antibody
<b>Family / Subfamily:</b>	Protein Kinase / PI3/PI4
<b>Host</b>	ATM antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	ATM antibody was raised against Human
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	ATM antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser1981 of human ATM.
<b>Specificity:</b>	Human ATM.
<b>Epitope:</b>	pSer1981
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Uses:</b>	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences