

ATM Rabbit anti-Human Polyclonal (pSer1981) (PE) Antibody - LS-C224865 - LSBio

CatalogID:	LS-C224865
Target:	ataxia telangiectasia mutated (ATM)
Synonyms:	ATM Antibody, A-T mutated Antibody, ATA Antibody, ATC Antibody, ATD Antibody, ATE Antibody, AT mutated Antibody, Ataxia telangiectasia Antibody, ATDC Antibody, Serine-protein kinase ATM Antibody, AT1 Antibody, Ataxia telangiectasia mutated Antibody, TEL1 Antibody, TELO1 Antibody
Family / Subfamily:	Protein Kinase / PI3/PI4
Host	ATM antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ATM antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	ATM antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser1981 of human ATM.
Specificity:	Human ATM.
Epitope:	pSer1981
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences