

DGKA Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C225113 - LSBio

CatalogID:	LS-C225113
Target:	diacylglycerol kinase, alpha 80kDa (DGKA)
Synonyms:	DGKA Antibody, 80 kDa diacylglycerol kinase Antibody, DAG kinase alpha Antibody, Diglyceride kinase alpha Antibody, DGK-alpha Antibody, Diacylglycerol kinase alpha Antibody, Diacylglycerol kinase, alpha Antibody, DAGK1 Antibody, DAGK Antibody
Family / Subfamily:	Non-protein Kinase / Diacylglycerol kinase
Host	DGKA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	DGKA antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	DGKA antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human DGKA.
Specificity:	Human DGKA
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences