

DGKQ Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C225126 - LSBio	
CatalogID:	LS-C225126
Target:	diacylglycerol kinase, theta 110kDa (DGKQ)
Synonyms:	DGKQ Antibody, DAGK7 Antibody, DAG kinase theta Antibody, DAGK4 Antibody, Diglyceride kinase theta Antibody, DAGK Antibody, DGK-theta Antibody, Diacylglycerol kinase theta Antibody
Family / Subfamily:	Non-Protein Kinase / Diacylglycerol kinase
Host	DGKQ antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	DGKQ antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	DGKQ antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human DGKQ.
Specificity:	Human DGKQ.
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	IHC, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	