

EPHA2 / EPH Receptor A2 Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C225345 - LSBio

CatalogID:	LS-C225345
Target:	EPH receptor A2 (EPHA2)
Synonyms:	EPHA2 Antibody, CTPA Antibody, EPH receptor A2 Antibody, Ephrin type-A receptor 2 Antibody, Epithelial cell kinase Antibody, ECK Antibody, Sek-2 Antibody, ARCC2 Antibody, Soluble EPHA2 variant 1 Antibody, CTPP1 Antibody
Family / Subfamily:	Protein Kinase / Ephrin Receptor
Host	EPHA2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	EPHA2 / EPH Receptor A2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	EPHA2 / EPH Receptor A2 antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human EphA2.
Specificity:	Human EPHA2.
Epitope:	N-Terminus
Reactivity:	Human, Mouse, Hamster
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences