

EPHA3 / EPH Receptor A3 Rabbit anti-Human Polyclonal (C-Terminus) (PE) Antibody - LS-C225389 - LSBio

CatalogID:	LS-C225389
Target:	EPH receptor A3 (EPHA3)
Synonyms:	EPHA3 Antibody, CEK4 kinase Antibody, EPH receptor A3 Antibody, Ephrin type-A receptor 3 Antibody, Eph-like tyrosine kinase 1 Antibody, ETK1 Antibody, HEK4 Antibody, EPH-like kinase 4 Antibody, HEK Antibody, Human embryo kinase 1 Antibody, MEK4 kinase Antibody, Human embryo kinase Antibody, Tyrosine-protein kinase TYRO4 Antibody, TYRO4 Antibody, EK4 Antibody, Eph kinase A3 Antibody, ETK Antibody
Family / Subfamily:	Protein Kinase / Ephrin Receptor
Host	EPHA3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	EPHA3 / EPH Receptor A3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	EPHA3 / EPH Receptor A3 antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human EphA3.
Specificity:	Human EPHA3
Epitope:	C-Terminus
Reactivity:	Human, Hamster
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences