

ICOS Rabbit anti-Human Polyclonal (C-Terminus) (PE) Antibody - LS-C225584 - LSBio

CatalogID:	LS-C225584
Target:	inducible T-cell co-stimulator (ICOS)
Synonyms:	ICOS Antibody, AILIM Antibody, CD278 Antibody, CVID1 Antibody, Inducible costimulator Antibody, Inducible T-cell co-stimulator Antibody, CD278 antigen Antibody, Inducible T-cell costimulator Antibody
Host	ICOS antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ICOS antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ICOS antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human ICOS.
Specificity:	Human ICOS.
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences