

IRAK1 / IRAK Mouse anti-Human Monoclonal (aa530-693) (FITC) Antibody - LS-C225690 - LSBio	
CatalogID:	LS-C225690
Target:	interleukin-1 receptor-associated kinase 1 (IRAK1)
Synonyms:	IRAK1 Antibody, IRAK-1 Antibody, Pelle Antibody, Il1rak Antibody, IRAK Antibody, Pelle homolog Antibody
Family / Subfamily:	Protein Kinase / IRAK(Pelle)
Host	IRAK1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IRAK1 / IRAK antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	IRAK1 / IRAK antibody was raised against recombinant protein corresponding to aa530-693 of human IRAK1.
Specificity:	Human IRAK1 / IRAK.
Epitope:	aa530-693
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	